






Note: Part numbers and standards' values may change. Consult LECO for the latest information.



See reverse side for contact information

**OPERATING SUPPLIES**

781-335 Quartz Boat  
782-059 Nickel Boat Liner (10/pkg)  
783-299 Combustion Tube

**CATALYSTS AND REAGENTS**

501-171-HAZ Anhydrone 454 g (1 lb.) bottle  
501-170 Copper Oxide 50 g bottle  
502-156 Fluorhib 454 g (1 lb.) bottle  
502-493 Fluorhib HP 35 g bottle  
501-081 Glass Wool 454 g (1 lb.) bottle  
502-174-HAZ Lecosorb 500 g bottle, (20–30 mesh)  
502-176-HAZ Lecosorb 500 g bottle, (8–20 mesh)  
501-608 Quartz Wool 200 g bottle

**REAGENT TUBES**

771-373 Glass Anhydrone Tube, requires two 772-520 O-rings  
782-428 Quartz Catalyst Heater Tube, requires two  
772-520 O-rings

**OPERATION ACCESSORIES**

780-044 Boat Puller Assembly  
761-929 Crucible Tongs  
601-483 Printer Paper  
601-484 Printer Ribbon  
782-508 Quartz Wool Puller/Packing Tool Assembly  
501-614 Spatula  
625-505-137 Quartz Boat Tweezer

**CALIBRATION SAMPLES**

501-034 12% Carbon (50 g); Calcium Carbonate  
502-029 1% Carbon (50 g); Synthetic Carbon  
502-030 5% Carbon (50 g); Synthetic Carbon  
502-091 12.3% Moisture (50 g); Calcium Oxalate

*NOTE: Part numbers and standards' values may change.  
Consult LECO for latest information.*

*Please Note: Quartz Tubes are typically utilized for higher temperature applications (i.e. catalyst tubes). Glass reagent/filter tubes are designed for lower temperature applications. Glass tubes should not be used in place of quartz tubes.*

**United States:**

Ph: 1-800-292-6141  
Fax: 269-982-8977

**Canada:**

Ph: 1-800-268-0215  
Fax: 905-564-6582

represented by:

**Other Countries:**

Ph: 269-983-5531  
Fax: 269-983-3850

© 2010 LECO Corporation